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QuickGEN First-Beer yeast and bacteria differentiation PCR Kit

Order No.: QTPYB 0096 white

(Version 02/18)

1. Intended use

Differentiation of beer spoilage bacteria and yeast in beer and beer mixing drinks.

The following bacteria and yeasts are differentiated:

Lactobacilli, *Pediococci*, Acetic acid bacteria, *Enterobacteriaceae*, Bottom-fermented yeast, Top-fermented yeast and wild yeast.

2. Test principle

The TaqMan[®] real-time PCR is based on hot-start PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ; ROX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

To avoid false negative PCR-results an Inhibition Control is amplified together in one tube with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **Tubes 4, 5, 7 and 8 contain lyticase.**

3. Kit contents

The kit contains sufficient reagents for 96 PCR reactions (24 samples):

1 x Premix QTPYB	white cap
12 x QTPYB mix (freeze-dried, tubes 4, 5, 7 und 8 contain lyticase)	tube strips
12 x cap strips	
1 x ddH ₂ O (for PCR negative controls)	colourless cap

4. Storage conditions

The QTPYB mix is freeze-dried in the strips and should be stored at 4 – 8 °C (39 – 46 °F).

Keep Premix at - 20 °C (- 4 °F). Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The strips with the QTPYB mix contain the fluorescent labelled probes and should be handled light protected.

All reagents are stable for 6 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR machine for low profile tubes white

Centrifuge for tube strips

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

optional: Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

Colour Compensation Kit (GEN-IAL® TPBDLC480 CC 0005)

LC480: Use a Colour compensation prior to first use of QTPYB kit

6. PCR

6.1. PCR-Setup

Before every use thoroughly mix Premix and centrifuge briefly. Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

PCR-reaction setup:

PCR-Components	Amount (µL)
Premix QTPYB	15.0
Sample-DNA	5.0*
Total volume	20.0

* if using the Simplex Easy DNA Kit: add 2.5 µL DNA and 2.5 µL sterile ddH₂O

Pipetting scheme:

Tube 1*	Premix 15.0 µL	ddH ₂ O 5.0 µL, Attention: pipette no sample-DNA (negative control)
Tube 2*	Premix 15.0 µL	ddH ₂ O 5.0 µL, Attention: pipette no sample-DNA (positive control)
Tube 3*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 4*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 5*	Premix 15.0 µL	sample-DNA of sample 1 5.0 µL
Tube 6*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL
Tube 7*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL
Tube 8*	Premix 15.0 µL	sample-DNA of sample 2 5.0 µL

* please notice the orientation of the tubes. They are marked on the strip at Tube 1 with an A and at Tube 8 with an H.

6.2 PCR-Program

6.2.1 PCR Program LC480

1. Click the button **tool** on the right side in the window **LightCycler 480 Software release 1.5.0. SP1**
2. Click the button **Detection formats** at the left side of the menu bar
3. Click **New** in the window **Detection formats** and name the experiment, for example: QTPYB
4. Open the window **Filter Combination Selection** and choose the following filter combinations: 465-510 / 533-580 / 533-610
5. Open the window **Selected Filter Combination** list and add the following amounts

Excitation Filter	Emission Filter	Name	Melt factor	Quant factor	Max. Integ. Time
465	510	465-510 FAM	1	10	2
533	580	533-580 HEX	1	10	2
533	610	533-610 610	1	10	2

6. **Close** the window
7. Click the button **New Experiment** on the right side of the menu bar
8. From the pull-down menu **Detection formats** choose the defined experiment, click the button **Customize** and check the detection formats. All of them have to be activated.
9. Choose **ok**
10. Write the following program:

1. Programm Name: **Heat**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
37	None	00:15:00	4.40		0	0	0
95	None	00:15:00	4.40		0	0	0

2. Programm Name: **Ampli**

Cycles **35** Analysis Mode **Quantification**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
95	None	00:00:15	4.40		0	0	0
62	Single	00:00:20	2.20		0	0	0

3. Programm Name: **Cool**

Cycles **1** Analysis Mode **None**

Target (°C)	Acquisition Mode	Hold (hh:mm:ss)	Ramp Rate (°C/s)	Acquisitions (per °C)	Sec Target (°C)	Step size (°C)	Step Delay (cycles)
40	None	00:00:20	2.20		0	0	0

Optional: Saving the program as **run template**. The template button allows to select and apply a template to the currently open object and to save the currently open object as a template. Click the clamp beside the button **Apply Template** . Click **Save as template** and save the file in **templates**.

11. Click the button **Subset editor** on the left side
12. Click the button **+** and **New Subset 1** appears
13. Mark the wells in the **New Subset 1 Settings** window
14. Click the button **Apply**
15. Click the button **Sample editor** on the left side
16. **Very important:** Activate In the window **Step 1 Select Workflow: Abs.Quant**
17. Choose the subset **New Subset** in the window **Step 2 Select Samples**
18. Define your probes in the **Sample table**
19. Switch to the button **Experiment** and start the run with the button **start run**

7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer. If the automatic data analysis of the completed run is not satisfying, the threshold can be adjusted manually above the background signals.

The negative control (NTC, tube 1) has to be negative in the FAM-, and HEX-channel and positive in the ROX-channel (Inhibition Control). Is the Inhibition Control in the ROX-channel negative and/or the results in the FAM-/ HEX-channel positive, PCR has to be repeated.

The Inhibition Controls in the samples (tubes 4/7) have to be positive in the ROX-channel. Is the Inhibition Control in the ROX-channel negative inhibitors or high amount of DNA are in the sample reaction.

The Positive Control (tube 2) has to be positive in the FAM-channel, if it is negative PCR, has to be repeated.

* Wild yeast 1 detects:

<i>Dekkera anomala</i>	<i>Dekkera bruxellensis</i>	<i>Dekkera custersiana</i>	<i>Dekkera naardenensis</i>	<i>Debaromyces hansenii</i>
<i>Hanseniaspora guilliermondii</i>	<i>Hanseniaspora osmophila</i>	<i>Hanseniaspora uvarum</i>	<i>Issotchenkia orientalis</i>	<i>Kazachstania exigua</i>
<i>Kluyveromyces marxianus</i>	<i>Metschnikowia pulcherrina</i>	<i>Pichia anomala</i>	<i>Pichia fermentans</i>	<i>Pichia membranaefaciens</i>
<i>Saccharomyces diastaticus</i>	<i>Saccharomyces ludwigii</i>	<i>Torulasporea delbrückii</i>		

* Wild yeast 2 detects:

<i>Candida glabrata</i>	<i>Candida albicans</i>	<i>Candida kefyr</i>	<i>Candida intermedia</i>	<i>Candida parapsilosis</i>
<i>Candida sake</i>	<i>Candida tropicalis</i>	<i>Naumovozyma dairenensis</i>	<i>Pichia guilliermondii</i>	<i>Zygosaccharomyces bailii</i>
<i>Zygosaccharomyces rouxii</i>				

PCR-Analysis:

please notice the orientation of the tubes. They are marked on the strip at Tube 1 with an A and at Tube 8 with an H.

Tube	Sample	FAM	HEX	ROX
1	NTC	-	-	Inhibition Control
2	PTC	Positive Control	-	-
3	1	Enterobacteriaceae	<i>Lactobacillus / Pediococcus</i>	<i>Pediococcus</i>
4	1	Wild yeast 1*	Bottom fermented yeast	Inhibition Control
5	1	Wild yeast 2*	Top fermented yeast	Acetic acid bacteria
6	2	Enterobacteriaceae	<i>Lactobacillus / Pediococcus</i>	<i>Pediococcus</i>
7	2	Wild yeast 1*	Bottom fermented yeast	Inhibition Control
8	2	Wild yeast 2*	Top fermented yeast	Acetic acid bacteria

NTC: negative Control PTC: positive Control

Example for interpretation:

If sample 1 shows a positive FAM-signal in tube 4, wild yeast 1 is identified. If the sample contains mixtures, you have positive signals in different channels.

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method. GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.