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(Version 02/18)

1. Intended use

Differentiation of wild yeast in beer and beer mixing drinks.

The following yeast are differentiated:

<i>Candida glabrata</i>	<i>Candida albicans</i>	<i>Candida kefir</i>	<i>Candida intermedia</i>	<i>Candida parapsilosis</i>
<i>Candida sake</i>	<i>Candida tropicalis</i>	<i>Debaromyces hansenii</i>	<i>Dekkera anomala</i>	<i>Dekkera bruxellensis</i>
<i>Dekkera custersiana</i>	<i>Dekkera naardenensis</i>	<i>Hanseniaspora guilliermondii</i>	<i>Hanseniaspora uvarum</i>	<i>Hanseniaspora osmophila</i>
<i>Kluyveromyces marxianus</i>	<i>Metschnikowia pulcherrima</i>	<i>Pichia anomala</i>	<i>Pichia membranaefaciens</i>	<i>Pichia fermentans</i>
<i>Pichia guilliermondii</i>	<i>Issotchenkia orientalis</i>	<i>Kazachstania exigua</i>	<i>Naumovozyma dairenensis</i>	<i>Saccharomyces diastaticus</i>
<i>Saccharomycodes ludwigii</i>	<i>Torulaspora delbrückii</i>	<i>Zygosaccharomyces rouxii</i>	<i>Zygosaccharomyces bailii</i>	

2. Test principle

The TaqMan[®] real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (FAM/DQ; HEX/DQ) which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument. To avoid false negative PCR-results an Inhibition control is amplified together in one tube with the specific sequence. The system contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit). **The tubes contain lyticase.**

3. Kit contents

The kit contains sufficient reagents for 96 PCR reactions (12 samples):

1 x Premix QTPWY	white cap
12 x QTPWY mix (freeze-dried, tubes contain lyticase)	tube strips
12 x cap strips	
1 x ddH ₂ O (for PCR negative controls)	colourless cap

4. Storage conditions

The QTPWY mix is freeze-dried in the strips and should be stored at 4 – 8 ° C (39 – 46 °F).

Keep Premix at - 20 °C (- 4 °F). Avoid loss of sensitivity by repeating freezing and thawing more than 3 times. For irregular use aliquot the Premix.

The tube strips with the QTPWY mix contain fluorescent labelled probes and should be handled light protected.

All reagents are stable for 6 months, if they are stored correctly.

5. Materials required, but not provided

5.1. Instruments

Real-time PCR machine for high profile tubes

Centrifuge for tube strips

Pipettes

“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O

sterile filter tips

optional: Uracil N-Glycosylase (0.01 U/μL added to the PCR reaction mix)

6. PCR

6.1. PCR-Setup

Before every use thoroughly mix Premix and centrifuge briefly. Remove the film from the needed tube strips and pipette the PCR-Components. After pipetting close the tube strips with the provided cap strips.

PCR-reaction setup:

PCR-Components	Volume (µL)
Premix QTPWY	17.5
Sample-DNA	2.5
Total volume	20.0

Pipetting schema:

Tube 1*	Premix 17.5 µL	ddH ₂ O 2.5 µL, Attention: pipette no sample-DNA (negative control)
Tube 2*	Premix 17.5 µL	ddH ₂ O 2,5 µL, Attention: pipette no sample-DNA (positive control)
Tube 3*	Premix 17.5 µL	DNA of sample 1 2.5 µL
Tube 4*	Premix 17.5 µL	DNA of sample 1 2.5 µL
Tube 5*	Premix 17.5 µL	DNA of sample 2 2.5 µL
Tube 6*	Premix 17.5 µL	DNA of sample 2 2.5 µL
Tube 7*	Premix 17.5 µL	DNA of sample 3 2.5 µL
Tube 8*	Premix 17.5 µL	DNA of sample 3 2.5 µL

* please notice the orientation of the tubes: they are marked with 1,2,3.....

6.2 PCR-Program

Note:

For the use of UNG the thermal cycler programm has to be changed according to manufacturers` instructions.

Step	Time	Temp.	
Lyticasetreatment	15 min	37 °C	
Initial denaturation of DNA	15 min	95 °C	
Cycling Denaturation	10 sec	95 °C	Cycle 40 x
Cycling Annealing/ Elongation	20 sec	60 °C	

7. Evaluation

The evaluation has to be made according to the data analysis program recommended by the real-time instrument manufacturer. If the automatic data analysis of the completed run is not satisfying, the threshold can be adjusted manually above the background signals.

The negative control (NTC, Tube 1) has to be negative in the FAM channel and positive in the HEX-channel (Inhibition Control). Is the Inhibition Control in the HEX-channel negative and/or the results in the FAM-channel positive, PCR has to be repeated.

The Inhibition Control of the sample has to be positive in the HEX-channel. Is the Inhibition-Control in the HEX-channel negative inhibitors or high amount of DNA are in the sample reaction.

The Positive Control (Tube 2) has to be positive in the FAM-channel. If it is negative, PCR has to be repeated.

PCR-Analysis:

please notice the orientation of the tubes: they are marked with 1,2,3.....

	Channels	Species	Signal
Tube 1 (NTC)	FAM	Negative control (NTC)	negative
	HEX	Inhibition-Control	positive
Tube 2 (PTC)	FAM	Positive control (PTC)	positive
	HEX	Inhibition-Control	
Tube 3 (Sample 1)	FAM	Wild yeast 1	
	HEX	Inhibition-Control	
Tube 4 (Sample 1)	FAM	Wild yeast 2	
	HEX	Inhibition control	
Tube 5 (Sample 2)	FAM	Wild yeast 1	
	HEX	Inhibition control	
Tube 6 (Sample 2)	FAM	Wild yeast 2	
	HEX	Inhibition control	
Tube 7 (Sample 3)	FAM	Wild yeast 1	positive
	HEX	Inhibition control	
Tube 8 (Sample 3)	FAM	Wild yeast 2	positive
	HEX	Inhibition control	

Example for interpretation:

Sample shows a positive FAM-signal in Tube 3: sample 1 is positive for wild yeast group 1

Wild yeast 1 detects:

<i>Dekkera anomala</i>	<i>Dekkera bruxellensis</i>	<i>Dekkera custersiana</i>	<i>Dekkera naardenensis</i>	<i>Debaromyces hansenii</i>
<i>Hanseniaspora guilliermondii</i>	<i>Hanseniaspora osmophila</i>	<i>Hanseniaspora uvarum</i>	<i>Issotchenkia orientalis</i>	<i>Kazachstania exigua</i>
<i>Kluyveromyces marxianus</i>	<i>Metschnikowia pulcherrima</i>	<i>Pichia anomala</i>	<i>Pichia fermentans</i>	<i>Pichia membranaefaciens</i>
<i>Saccharomyces diastaticus</i>	<i>Saccharomycodes ludwigii</i>	<i>Torulaspora delbrückii</i>		

Wild yeast 2 detects:

<i>Candida glabrata</i>	<i>Candida albicans</i>	<i>Candida kefyr</i>	<i>Candida intermedia</i>	<i>Candida parapsilosis</i>
<i>Candida sake</i>	<i>Candida tropicalis</i>	<i>Naumovozyma dairenensis</i>	<i>Pichia guilliermondii</i>	<i>Zygosaccharomyces bailii</i>
<i>Zygosaccharomyces rouxii</i>				

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.