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2. Experiment: Initial denaturation
Analysis Mode: None

Cycles: 1

Target Temperature (°C)	Incubation Time (hrs:min:sec)	Temperature Transition Rate (°C/s)	Acquisition Mode
95	00:10:00	20.00	None

3. Experiment: Amplification

Analysis Mode: Quantification

Cycles: 35 - 45

Target Temperature (°C)	Incubation Time (hrs:min:sec)	Temperature Transition Rate (°C/s)	Acquisition Mode
94	00:00:05	20.00	None
60	00:00:25	20.00	Single

4. Experiment: Cool
Analysis Mode: None

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
Art. No.: TPYPM 0050

1. Intended use

Detection of *Pichia membranaefaciens* in beverages by real-time PCR.

2. Test principle

The real-time PCR is based on hot-start-PCR and sequence-specific dual labelled probes (TaqMan®), which, when accurately hybridised, emit a measurable fluorescent signal of a defined wavelength in the extension phase. The increase of signal is continuously measured in a real-time PCR detection instrument.

The kit contains one specific system for the detection of *Pichia membranaefaciens*. The system emits a maximum fluorescent signal at 518 nm (FAM-channel) and contains dUTP. Optional: the use of Uracil-N-Glycosylase will eliminate any contamination with Uracil containing amplicons from former PCRs (the enzyme is not part of this kit).

3. Kit contents

The kit contains sufficient reagents for 50 PCR reactions.

1 x Premix TPYPM	white cap
1 x TPYPM mix (freeze-dried)	dark vial, red cap
1 x ddH ₂ O (for solving dried reagents)	colourless cap
1 x Control-DNA (<i>P. membranaefaciens</i> , freeze-dried)	yellow cap

4. Storage conditions

The TPYPM mix and the Control-DNA are freeze-dried, they have to be solved in ddH₂O prior to use (see 6.1).

All freeze-dried reagents should be stored cool at 4 – 8 °C (35 – 46 °F), the Premix at – 20 °C (– 4 °F).

Once the freeze-dried reagents are dissolved they are stable for 4 weeks at 4 – 8°C (35 – 46 °F) or up to 6 months at –20 °C (–4 °F). Prevent freeze/thaw cycles. Prepare aliquots, if the PCR kit should be used for a longer period than 4 weeks. Frozen aliquots (–20°C) (–4°F) are stable for 6 months.

The TPYPM mix (red cap) contains the fluorescent labelled probe and should be handled light protected.

5. Materials required, but not provided

5.1. Instruments

Realtime PCR instrument with at least FAM-channel
Standard benchtop mini-centrifuge (1.5 – 2.0 mL vials)
Pipettes for 0.5 – 1000 µL
“Vortex”

5.2. Reagents and plastic ware

sterile ddH₂O
sterile reaction vessels 1.5 – 2.0 mL
sterile optical tubes (plastic) or glass capillaries (LightCycler™)
sterile filter tips
optional: Uracil N-Glycosylase (0.01 U/µL added to the PCR reaction mix)

6. PCR

6.1. PCR-Setup

Preparations before first use:

If you use the kit for the first time, the freeze-dried kit components have to be shortly centrifuged and carefully resolved:

- add 80 µL sterile ddH₂O to freeze-dried TPYPM mix
- add 55 µL sterile ddH₂O to the freeze-dried Control-DNA
- after 15 minutes mix well

Before every use thoroughly mix all PCR components and centrifuge briefly.

PCR-reaction per sample:

PCR-components	amount (µL)
Premix TPYPM	13.5
TPYPM mix (dissolved in ddH ₂ O)	1.5
Sample or Control-DNA	up to 5.0
ddH ₂ O	ad. 20
total volume	20.0

Prepare a master mix by mixing Premix and TPYPM mix

1. Multiply said volumes with the number of PCR preparations including controls (positive control, negative control, extraction control), taking into account pipette reserves of approximately 5-10 %.
2. Divide 15 µL of the PCR-master mix among the individual reaction vessels, making sure that, prior to the first filling, the tip of the pipette has been moistened.
3. Add the sample DNA using a fresh tip with each DNA filling and mix it with the master mix at the same time, (with DNA volume less than 5 µL add ddH₂O to the end volume of 20 µL per reaction). Add 5 µL of the Control-DNA and 5 µL of ddH₂O for the negative control reaction.
4. Close the tubes immediately and centrifuge them according to the recommendations of the real-time PCR instrument manual.
5. Place tubes/plate/capillaries into the thermo block/carousel of the instrument, close lid and start run.

Work swiftly to avoid warming up and keep away from light.

6.2 PCR-Programme

Note:

For the use of UNG the thermal cycler programmes have to be changed according to manufacturers' instructions.

6.2.1 Real-time instruments using optical tubes or 96-well plates

The Premix contains no dye for normalising the fluorescence signal

Step	Time	Temp.	Acquisition Mode	Status
Initial denaturation of DNA	10 min	95 °C	none	Hold
Cycling Denaturation	15 sec	95 °C	none	Cycle 35 - 40
Cycling Annealing/ Elongation	30 sec	60 °C	single	
Cooling	15 sec	40 °C	none	Hold

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7. Interpretation of results

The evaluation has to be made according to the data analysis programme recommended by the real-time instrument manufacturer. If the automatic data analysis of the completed run is not satisfying, the threshold can be adjusted manually above the background signals.


The *Pichia membranaefaciens* specific signal is measured at:

FAM calibrated channel at 518 nm (in a Rotor-Gene™ channel 520 nm, in a LightCycler™ 2.0 channel 530 nm)

A sample is positive, if there is a detectable fluorescence increase in the FAM/530 nm channel and the negative controls show no amplification. The positive controls should have a positive fluorescence signal.

A sample is negative, if there is no detectable fluorescence increase in the FAM/530 nm channel and the positive controls have a positive fluorescence signal. The negative controls show no amplification in FAM e.g. 530 channel.

analysis flowchart

Note:

The polymerase-chain reaction (PCR) is protected by patents and requires a licence from Hoffmann-LaRoche Inc.. The provided product does not authorise the purchaser for the commercial use of this method.

GEN-IAL makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, GEN-IAL will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. GEN-IAL shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

For research purposes only.