RIDASCREEN[®] Sulfamethazin

Art. No. R3011

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RIDASCREEN[®] Sulfamethazin

Brief information

RIDASCREEN[®] Sulfamethazin (Art. No.: R3011) is a competitive enzyme immunoassay for the quantitative determination of sulfamethazine in milk, meat, honey, liver, kidney, fish, shrimp and egg.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	<u>milk</u> : direct use <u>meat</u> : homogenization, centrifugation <u>honey</u> : extraction, clean-up with RIDA [®] C18 column, evaporation, reconstitution, precipitation <u>liver, kidney, fish</u> : homogenization, centrifugation, precipitation <u>shrimp, egg</u> : homogenization, extraction, evaporation, reconstitution
Time requirement:	sample preparation (for 10 samples) according to sample matrix..... 5 min to 2 h test implementation (incubation time) 45 min
Limit of Detection (corresponding to the standard substance)	milk..... approx. 4 µg/L meat (bovine, porcine)..... approx. 5 µg/kg meat (poultry) approx. 10 µg/kg honey* approx. 10 µg/kg liver approx. 6 µg/kg kidney..... approx. 10 µg/kg fish approx. 7 µg/kg shrimp approx. 15 µg/kg egg approx. 16 µg/kg

* For honey the LOD can deviate due to matrix effects.

Recovery	milk.....	approx. 110 %
(corresponding to the	meat (bovine, porcine).....	approx. 94 %
standard substance)	meat (poultry)	approx. 115 %
	honey	approx. 88 %
	liver	approx. 89 %
	kidney.....	approx. 96 %
	fish	approx. 88 %
	shrimp	approx. 111 %
	egg.....	approx. 107 %

The specificity of the RIDASCREEN® Sulfamethazin test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from the values determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the limit of detection and the recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Specificity:

in buffer system	Sulfmethazine (standard substance)	100 %
	Sulfamerazine	approx. 4 %
	Sulfacetamide, sulfachloropyridazine, sulfachloropyridazine -(phenyl-13C6), sulfadiazine, sulfadimethoxin, sulfadoxine, sulfaguanidine, sulfamethoxazole, sulfamethoxypyridazine, sulfamethizole, sulfamoxole, sulfanilamide, sulphaphenazole, sulfapyridine, sulfaquinoxaline, sulfisoxazole.....	< 1%

In order to increase the quality of assessment when performing ELISA-procedures, we refer additionally to our Good ELISA Practice (GEP) – Manual in the respective version. Such lists minimum standards concerning the framework conditions when using test-kits of R-Biopharm AG and performing ELISA-analysis. The manual can be retrieved, printed and downloaded under the website www.r-biopharm.com/products/food-feed-analysis.

Related products

RIDA [®] Sulfamethazin Spiking Solution	(R3098)
RIDASCREEN [®] Sulfonamide	(R3004)
RIDA [®] Sulfonamide-Sulfamethoxypyridazin spiking solution	(R3099)

1. Intended use

RIDASCREEN[®] Sulfamethazin is a competitive enzyme immunoassay for the quantitative determination of sulfamethazine in milk, meat, honey, liver, kidney, fish, shrimp and egg.

2. General

Sulfonamides are widely used as feed additives, mainly for fattening of calves and pigs. Combined with inhibitors of dihydrofolate reductase such as trimethoprim, tetroxoprim, or pyrimethamine sulfonamides are also used in veterinary medicine for the treatment of intestinal infections, mastitis, pneumonitis and other (systemic) diseases. Sulfonamide residues may therefore occur in food of animal origin such as meat and milk. Particularly the transmission of the carcinogenic sulfamethazine represents a threat to human health. The EU established a maximum residue limit (MRL) for all substances of the sulfonamide-group of 100 µg/kg (ppb) in muscle, fat, liver, kidney and milk.

3. Test principle

The basis of the test is the antigen-antibody reaction. The microtiter wells are coated with capture antibodies directed against anti-sulfamethazine antibodies. The microtiter plate is already functionalized with anti-sulfamethazine antibodies, which are bound by the immobilized capture antibodies. Standards or sample and sulfamethazine conjugate are added into the wells. Free and enzyme conjugated sulfamethazine competes for the antibody binding sites (competitive enzyme immunoassay). Any unbound conjugate is then removed in a washing step. The conjugate converts the chromogen into a blue product. The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the sulfamethazine concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 Wells
Standard 1	white	Ready to use	0 µg/L	1.3 ml
Standard 2	white	Ready to use	3 µg/L	1.3 ml
Standard 3	white	Ready to use	9 µg/L	1.3 ml
Standard 4	white	Ready to use	27 µg/L	1.3 ml
Standard 5	white	Ready to use	81 µg/L	1.3 ml
Standard 6	white	Ready to use	162 µg/L	1.3 ml
Wash buffer salt Tween		Salt for dissolving		
Conjugate	red	Ready to use		7.5 ml
Substrate/Chromogen Red Chromogen Pro	brown	Ready to use		10 ml
Stop solution	yellow	Ready to use		14 ml

5. Materials required but not provided

5.1. Equipment:

Equipment	milk	meat	honey	liver	kidney	fish	shrimp egg
microtiter plate spectrophotometer (450 nm)	•	•	•	•	•	•	•
graduated pipettes		•	•	•	•	•	•
variable 5 - 10 µl micropipettes			•		•	•	
variable 20 - 200 µl and 200 - 1000 µl micropipettes	•	•	•	•	•	•	•
mixer		•		•	•	•	•
shaker		•		•	•	•	•
centrifuge		•	•	•	•	•	•
evaporator			•				•
vortex	•	•	•	•	•	•	•
RIDA® C-18 column			•				

5.2. Reagents:

Reagent	milk	meat	honey	liver	kidney	fish	shrimp egg
extraction buffer			•				
methanol			•				•
Carrez I			•	•	•	•	
Carrez II			•	•	•	•	

Extraction buffer for honey: 50 mM sodium acetate buffer, pH 5

- Dissolve 4.1 g sodium acetate (CH_3COONa) in 900 ml distilled water, adjust pH to 5 with HCl and fill up to 1000 ml with distilled water

Carrez I:

- Carrez I: 1.52 g Potassium ferrocyanide(II) $\times 3 \text{ H}_2\text{O}$; fill up to 10 ml with distilled water

Carrez II:

- Carrez II: 2.99 g Zinc sulfate $\times 7 \text{ H}_2\text{O}$; fill up to 10 ml with distilled water

6. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C. Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C.

The substrate/chromogen is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after expiry of the kit (see kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of deterioration of reagents

- any bluish coloration of the reddish substrate/chromogen prior to test implementation
- a value of less than 0.6 absorbance units ($A_{450\text{ nm}} < 0.6$) for the zero standard

9. Preparation of Samples

The samples should be stored in a cool place, protected against light.

9.1. Milk

- vortex a representative sample amount for homogenization
- use 50 µl milk per well in the assay

9.2. Meat

- remove fat from a representative amount of meat and homogenize the meat with an Ultra-Turrax or a mixer
- add 2 ml washing buffer (see 10.1.) to 1 g of homogenized sample
- vortex for 10 seconds and then shake vigorously upside down for 5 minutes
- centrifuge: 10 min / 4000 g / room temperature (20 - 25 °C)
- if there is an upper fat layer, remove it completely (e.g. with a pasteur pipette)
- use 50 µl of the supernatant per well in the assay

9.3. Honey

- dissolve 1 g of honey in 6 ml of extraction buffer (see 5.2.)
- centrifuge: 10 min / 3,000 g / room temperature (20 - 25 °C)

Purify all of the sample extract with RIDA[®] C-18 column (Art. No. R2002):

- wash the column with 2 ml methanol
- equilibrate the column with 2 ml of extraction buffer
- apply all of the sample extract and wash with 8 ml extraction buffer
- completely remove any fluid from the column by flushing the column with N₂ or air
- elute sample slowly with 1 ml of 100 % methanol (approx. 15 drops per second)
- evaporate the eluate to dryness under an nitrogen or air stream at 60 °C
- reconstitute the dry residue in 1 ml washing buffer
- add 10 µl of Carrez I, vortex
- add 10 µl of Carrez II, vortex
- centrifuge: 10 min / 4,000 g / room temperature
- use 50 µl of the supernatant per well in the assay

9.4. Liver

- remove fat from a representative amount of sample and homogenize the liver with an Ultra-Turrax or a mixer
- add 2 ml washing buffer (see 10.1.) to 1 g of homogenized sample
- vortex vigorously for 10 seconds and then shake upside down for 5 minutes
- centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- to 1 ml of the supernatant add 20 µl Carrez I, vortex
- add 20 µl Carrez II, vortex
- centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- use 50 µl of the supernatant per well in the assay

9.5. Kidney

- remove fat from the sample and homogenize the kidney with an Ultra-Turrax or a mixer
- add 2 ml washing buffer (see 10.1.) to 1 g of homogenized sample
- vortex for 10 seconds and then shake vigorously upside down for 5 minutes
- centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- to 1 ml of the supernatant add 10 µl Carrez I, vortex
- add 10 µl Carrez II, vortex
- centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- use 50 µl of the supernatant per well in the assay

9.6. Fish

- remove fat and skin from a representative amount of the sample and homogenize the fish with an Ultra-Turrax or a mixer
- add 2 ml washing buffer (see 5.1.) to 1 g of homogenized sample
- vortex for 10 seconds and then shake vigorously upside down for 5 minutes
- centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- to 1 ml of the supernatant add 5 µl Carrez I, vortex
- add 5 µl Carrez II, vortex
- centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- use 50 µl of the supernatant per well in the assay

9.7. Shrimp/egg

- homogenize the sample with an Ultra-Turrax or a mixer
- add 2 ml methanol to 1 g of homogenized sample
- vortex vigorously for 30 seconds
- centrifuge: 10 min / 4,000 g / room temperature (20 - 25 °C)
- transfer 1 ml of the supernatant into a new vial and evaporate to complete dryness under nitrogen or air stream at 60 °C

- reconstitute the dried residue in 1 ml washing buffer (see 10.1.)
- use 50 µl of the supernatant per well in the assay

10. Test implementation

10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

As **wash buffer** a PBS tween buffer is needed. Please use the wash buffer salt (see 4.) contained in the kit. Dissolve the entire buffer salt in one liter of distilled water. The ready to use washing buffer expires after approx. 4 - 6 weeks at 2 - 8 °C (36 - 46 °F).

Alternative: Dissolve the contents of the pouch in 100 ml of distilled water to obtain a 10fold concentrated washing buffer. This solution expires after approx. 8 - 12 weeks, store at room temperature (20 - 25 °C). Use 1 part 10fold concentrate and dissolve with 9 parts of distilled water to obtain the ready to use washing buffer.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run in duplicate. Record standard and sample positions.
2. Add 50 µl of each standard or prepared sample to separate duplicate wells.
3. Add 50 µl of the conjugate to the bottom of each well, mix gently by shaking the plate manually and incubate for 30 min at room temperature (20 - 25 °C).
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
5. Add 100 µl of substrate/chromogen solution to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C) in the dark.
6. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 30 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win.NET (Art. No. Z9996), is available to evaluate the RIDASCREEN[®] enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Remark for the calculation without software:

$$\frac{\text{absorbance standard (or sample)}}{\text{absorbance zero standard}} \times 100 = \% \text{ absorbance}$$

The zero standard is thus made equal to 100 % and the absorbance values are quoted in percentages. The values calculated for the standards are entered in a system of coordinates on semilogarithmic graph paper against the sulfamethazine concentration [µg/L].

In order to obtain the sulfamethazine concentration in µg/L / µg/kg (ppb) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the regulation stated, the dilution factors are as follows:

milk	1
meat.....	3
honey	1
liver	3
kidney	3
fish	3
shrimp	3
egg.....	3

For further information or applications please contact your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG